


Distributed by
CARL REINER ■


Catheter Solutions for Jet Ventilation



Latex-Free
Phthalate-Free
PVC-Free


Jet Catheters, 1 Lumen, not Laser Resistant

PEDIATRIC REF CTFS-513-000
ADULT REF CTFS-517-000



Jet Catheters, 4 Lumen, not Laser Resistant

ADULT REF CTFS-546-010



NOT LASER RESISTANT

Lumen	Outer Diameter	Length	Patient Gas Analysis	Pressure Measurement	Package Contains	REF	Adapter REF
1 Lumen	2.7 mm	35 cm	no	PP	100 pieces	CTFS-513-000	CTFS-500-400
1 Lumen	4 mm	45 cm	no	PP	100 pieces	CTFS-514-000	CTFS-500-400
1 Lumen	4 mm	70 cm	no	PP	100 pieces	CTFS-517-000	CTFS-500-400
3 Lumen	5.8 mm	45 cm	no	Paw	5 pieces	CTFS-534-000	CTNI-610-400
4 Lumen	6.6 mm	45 cm	yes	Paw	5 pieces	CTFS-544-010	CTNI-610-400
4 Lumen	6.6 mm	65 cm	yes	Paw	5 pieces	CTFS-546-010	CTNI-610-400

LASER RESISTANT

Lumen	Outer Diameter	Length	Patient Gas Analysis	Pressure Measurement	Package Contains	REF	Adapter REF
1 Lumen	3.8 mm	45 cm	no	PP	5 pieces	CTFS-514-001	CTFS-500-400
1 Lumen	3.8 mm	65 cm	no	PP	5 pieces	CTFS-516-001	CTFS-500-400
2 Lumen	3.8 mm	40 cm	yes	PP	5 pieces	CTFS-524-011	CTNI-610-400
2 Lumen	3.8 mm	65 cm	yes	PP	5 pieces	CTFS-526-011	CTNI-610-400

For more details feel free to contact us

Manufactured by Medicoplast, Saarland